




# Useddairyequipment


## Laboratory homogeniser / Labor-Homogenisator / Homogénéisateur de laboratoire / Laboratorium homogenisator / Homogeneizador de laboratorio / Лабораторный гомогенизатор


Offer number	: 24-31-03
Make	: Bos
Type	: Laboratory homogeniser
Power	: 1,1 kW - 3,2 A - 680 rpm - 50 Hz.
Dimensions ( l x w x h )	: 52 x 50 x 67 cm.
Weight	: 118 kg's
Condition of the machine	: Good
Location	: The Netherlands
Terms of delivery	: FCA





 Laboratory homogenizer. 2 stages. Manual pressure control. 3 Pistons. Stainless steel clad. Current motor power: 1,1 kW - 3,2 A - 680 rpm - 50 Hz. Manometer on first stage. Total dimensions: 52 x 50 x 67 cm. Weight: 118 kg's

 Labor-Homogenisator. 2-stufig. Manuelle Druckregelung. 3 Kolben. Edelstahl verkleidet.. Aktuelle Motorleistung: 1,1 kW - 3,2 A - 680 U/min - 50 Hz. Manometer an der ersten Stufe. Gesamtmaße: 52 x 50 x 67 cm. Gewicht: 118 kg

 Homogénéisateur de laboratoire. 2 étages avec contrôle manuel de la pression. 3 pistons. Revêtement en acier inoxydable. Puissance moteur actuelle : 1,1 kW - 3,2 A - 680 tr/min - 50 Hz. Manomètre sur premier étage. Dimensions totales : 52 x 50 x 67 cm. Poids : 118 kg

 Laboratorium homogenisator. 2-traps. Handmatige drukregeling. 3 zuigers. RVS behuizing.. Huidig motorvermogen: 1,1 kW - 3,2 A - 680 tpm - 50 Hz. Manometer op de eerste trap. Totale afmetingen: 52 x 50 x 67 cm. Gewicht: 118 kg

 Homogeneizador de laboratorio. 2 etapas. Control de presión manual. 3 pistones. Revestido de acero inoxidable. Potencia actual del motor: 1,1 kW - 3,2 A - 680 rpm - 50 Hz. Manómetro en primera etapa. Dimensiones totales: 52 x 50 x 67 cm. Peso: 118 kg

 Лабораторный гомогенизатор. 2 ступени с ручным регулированием давления. 3 поршня. Покрытие из нержавеющей стали. Текущая мощность двигателя: 1,1 кВт - 3,2 А - 680 об/мин - 50 Гц. Манометр на первой ступени. Общие размеры: 52 x 50 x 67 см. Вес: 118 кг

**For more information check our website: [www.useddairyequipment.com](http://www.useddairyequipment.com)**

PLEASE NOTE: We do not intentionally send our offers to any unwilling recipients.

If you have been incorrectly included in our mailing list, or do not wish to receive further offers from us, please reply with **REMOVE** as subject, and use the same mail address that we have send to. Please accept our apologies for any inconvenience caused!